


ZAHNER

PRECISION IN ELECTROCHEMISTRY



ZENNIUM XC - Made in Germany

ZENNIUM XC


THE COMPACT POTENTIOSTAT

High
Dynamic
Precision


Technology

Main Specifications

- EIS frequency range 10 μHz – 5 MHz
- 32-bit DC and 24-bit AC resolution
- $\pm 5\text{ V}$ / $\pm 14\text{ V}$ voltage range
- $\pm 2\text{ A}$ over 12 current ranges
- Online data processing for outstanding EIS



Slow CV scans with a scan rate of 10 $\mu\text{V/s}$ on a highly capacitive system with the ZENNIUM potentiostat (32-bit DAC resolution) and a conventional potentiostat (16-bit DAC resolution).




Single frequency (1 Hz), single period impedance measurements on a 25 $\mu\Omega$ resistor vs. time. The measurement is carried out with 1 A amplitude.

” THE HIGH-END POTENTIOSTAT “


Our Strengths

Zahner potentiostats can carry out advance electrochemical measurements like NFRA and intelligent multi sine EIS measurement besides traditional electrochemical measurements.

For more information:




NFRA: ZENNIUM potentiostats can measure harmonics during the impedance measurement, making it capable to carry out non-linear frequency response analysis.



IM-Sine: ZENNIUM potentiostats can carry out intelligent multi-sine EIS measurements, significantly decreasing the total measurement time.

Accuracy Contour Plot

- $Z > 0.1\ \Omega$: potentiostatic mode, amplitude 10 mV
- $Z > 1\ \text{M}\Omega$: potentiostatic mode, amplitude 50 mV, shielded
- $Z < 0.1\ \Omega$: galvanostatic mode, amplitude 100 mA
- $Z < 0.01\ \Omega$: galvanostatic mode, amplitude 1 A
- Without DC bias voltage/current
- Specified at the BNC terminals



Specifications

Potentiostatic modes	potentiostatic, galvanostatic, pseudo-galvanostatic, rest potential, ZRA, off
ADC resolution	32 bit
Function generator	digital (analog: option ADF for scan rates up to 10 kV/s)
Harmonic reject	> 60 dB @ ½ full scale
Cell connection	2-, 3-, 4-terminal kelvin
Ground reference	grounded, floating

Frequency generator & analyzer	Low range	High range
EIS frequency range	10 µHz to 5 MHz	
AC amplitude	0 to 2 V, 24 bit resolution	0 to 6 V, 24 bit resolution
Accuracy	< 0.0025%	
Resolution	0.0025%, 10,000 steps/decade	

Output potentiostatic	Low range	High range
Controlled voltage	±5 V	±14 V
Resolution	2.5 nV	7.5 nV
Accuracy	±200 µV ± 10 ppm of reading	±600 µV ± 10 ppm of reading
Integral nonlinearity	typ. 4 ppm, max. 8 ppm	typ. 12 ppm, max. 24 ppm
Compliance voltage	±14 V	±14 V
Bandwidth	DC to 6 MHz @ 33 Ω load	
IR compensation	auto AC impedance technique, range 0 to 10 MΩ, resolution 0.012%	
Small signal rise time	150 ns to 200 µs in 5 steps, automatic selection	
Slew rate	15 MV/s	
Phase shift	10° @ 500 kHz	

Output galvanostatic	
Controlled current	±2 A
Current range	±1.9 nA to ±2 A in 12 current ranges
Resolution	32 bit ± 0.2 ppb of FS
Accuracy	±0.1% of reading ± 0.04% of FS, ≥1 µA to 100 mA ±0.4% of reading ± 0.2% of FS, < 1 µA or > 100 mA

Input	Low range	High range
Max. Input voltage	±5 V	±14 V
Voltage resolution	2.5 nV	7.5 nV
Voltage accuracy	±100 µV ± 5 ppm of reading	±300 µV ± 10 ppm of reading
DC current resolution	2 nA (32 bit)	
DC current accuracy	±0.05% of reading ± 0.04% of FS @ 1 µA ... 100 mA ±0.5% of reading ± 0.4% of FS @ 100 mA ... 2 A ±0.5% of reading ± 0.4% of FS @ 10 nA ... 1 µA ±0.5% of reading ± 125 fA @ < 1 nA (HiZ-Probe)	
Input impedance	> 10 TΩ ±5 pF typ. (Main) / > 1000 TΩ ±1 pF typ. (HiZ-Probe)	
Input leakage current	< ±200 fA typ., < ±5 pA max., / < ±10 fA typ. (HiZ-Probe)	
Impedance range	1 mΩ to 10 GΩ / 2% (Main) 100 mΩ to 10 MΩ / 0.2% 100 mΩ to 100 GΩ / 2% (HiZ-Probe) 10 µΩ to 1 GΩ / 2% (Gal) 1 mΩ to 10 MΩ / 0.2%	
Common mode rejection	> 86 dB @ 10 µHz to 100 kHz > 66 dB @ 100 kHz to 5 MHz	
Input channel phase-tracking acc.	±0.05° @ 10 µHz to 100 kHz ±0.125° @ 100 kHz to 5 MHz	
Equivalent effective input noise	1 µV rms / 100 fA rms @ 1 mHz to 10 Hz	



Remote integration possible via Python and C++. Check out complete API documentation.

PC interface	USB 2.0
Dimensions / Weight	160 × 255 × 385 mm ³ / 8 kg
Power supply	100/115/230 VAC, 50/60 Hz
Ambient temperature / humidity	+10 °C to +30 °C / < 60% without derating

Zahner Analysis


EIS fitting

- Create equivalent electrical circuits
- Fit impedance spectra
 - > Single fit
 - > Series fit
- ZHIT tool
- Significance plot
- Fitting accessible via API


Other techniques

- Cyclic voltammetry
 - > Peak determination
 - > Charge integration
- Tafel slope measurements
- Butler-Volmer measurements
- Analysis of photoelectrochemical measurements


Check out
Zahner
Analysis
videos:




Create your own equivalent electric circuit for EIS fitting



Impedance spectrum (Bode plot) of a battery with the equivalent electrical circuit



Impedance spectrum (Nyquist plot) with two time constants



CV scans measured at different scan speeds

ZHIT

The Zahner Analysis software features the unique ZHIT tool, which helps identifying artifacts in impedance spectra and allows reconstruction of artifact-free impedance spectra for fitting.

Significance Plot

Zahner Analysis software features an exclusive tool called the significance plot, which evaluates the frequency-dependent significance of equivalent circuit elements in the fitting.

Zahner-Elektrik
GmbH & Co. KG
contact@zahner.de
www.zahner.de

